

































































Compatibilidade	Conector	Modelo	Nº de vias
Abbott Medical	 AAMI (6 pin - Macho)	AAEM; AAXM; ABBV; ABEM; ABTM; ABBY Séries; ARTM; DHS; EK40; EK50; EKP-1; EKR-50; EO-8; EX-50.	3 / 5
AD Instruments	 AAMI (6 pin - Macho)	MAC LAB (ECG).	3
Air Shields / Isolette	 AAMI (6 pin - Macho)	HRRM-71; HRRM72; System; Neo Natal; System 5 New.	3 / 5
Anamed	 AAMI (6 pin - Macho)	AM-78100A; AM-78100B; AM-78600; Vital Lite; Vital Wave; AM-7810A; HP- 14056B.	3 / 5
	 DIN 60°	ECG AM – 1.	5
Artema / S&W	 P 10	CB-600.	3
Bese	 SV800	Biomonitor 1.	3 / 5
	 AMP (6 pin - Fêmea)	Biomonitor 4.	3
	 SV800	Biomonitor 6.	3 / 5
	 AMP (6 pin - Fêmea)	Biomonitor 7.	3 / 5
Berger	 SV500	CD-181; CD-188.	5
Biosys	 Redel 6pin 1 guia	Guardian; BPM300; BPM200; BPM808; Delta Plus.	3 / 5
Bircher	 DB15 (15 pin-macho)	365; 375; Cardiotracer; 375-11.	10
Burdick / Spacelabs	 DB15 (15 pin-macho)	Eclipse 4; Eclipse 400; Eclipse 8; Eclipse 800; Eclipse 850; Eclipse 850i; Eclipse Plus; Eclipse LE II BURDICK (com n/s começando em: 1067 ou 92304).	10
	 AAMI (6 pin - Macho)	900 Séries; Medic-3; Medic-4; Medic-5.	3 / 5

	 DB15 (15 pin-macho)	E350; E350i; E550; E560; EK10; Eclipse LE/LEI; Elite; Elite II.	10
Cambridge	 DB15 (15 pin-macho)	SC1000; MC6000; VS 6.	10
Cardioline	 DB15 (15 pin-macho)	E1; 1 R/3; DELTA 3/300/360.	10
	 SV500	E – 1.	5
Cmos Drake	 SV500	ISO-212.	5
		MDK-2000; MDK-2001.	4
	 DB9 - (Fêmea)	Life 400 Plus.	5
	 AAMI (6 pin - Macho)	2010.	5
	 SV600	MDK - 2007/2010.	5
Criticare System International (CSI)	 AAMI (6 pin - Macho)	506; 507; 1100; 1508; 507EL; 507ELC; 508; 508E; 508EP; 508ER; 8100.	3 / 5
		504US; 602-4USP; Scholar.	3 / 5
Critikon	 AAMI (6 pin - Macho)	Dinamap Plus; 8720; 9720.	3 / 5
Datascope	 AAMI (6 pin - Macho)	2000; 2001; 2002; 2100; 2200; 3000; 3000A; Balloon Pump; CS100; M/D3; M/D3A; M/D4; Passport (todos); Spectrum; System (Todos); Trio; 851; 861; 870; 870A; 871; 990; M/D2; M/D2J; M/D2S.	3 / 5
Datex / Ohmeda	 Datex 10pin Fêmea	Cadiocap I; Cadiocap II; Cadiocap 5; AS/3; CS/3; S/5; Sattelite Plus; IMM Light.	3 / 5
Del mar	 DB15 (15 pin-macho)	910; 920; 930; 940; 942; 950; 960; 970; 980.	10
Digicare		Digimax 5000.	3
	 SV500	Digimax 5500.	5
		Cardiotrace; Cardiotrace; Vitaltrace.	3
Dixtal	 AAMI (6 pin - Macho)	DX2010; DX2021.	3 / 5

	 DB9 - (Fêmea)	SDM-2000.	5
	 DB9 - (Fêmea)	DX-910/920.	3 / 5
	 DB15 (15 pin- fêmea)	Eletropágina: EP-3.	10
Drager	 AAMI (6 pin - Macho)	Infinity Vista XL; VitalERT 2000; 1000.	3 / 5
Emai / Transmai	 SV500	RX-10; RX-20.	3 / 5
	 AAMI (6 pin - Macho)	RX-300; MX-10; MX-20; MX-300.	3 / 5
Ergoline	 DB15 (15 pin-macho)	Autoruler, Autoscript 6/12; Cardiette, EK 3003 / 3012, Ergoline	10
Ecafix	 SV500	ECG 5; ECG 6.	5
	 AAMI (6 pin - Macho)	Maximo	3 / 5
	 AMP (6 pin - Fêmea)		3 / 5
	 SV500		3 / 5
	 SV500	TC-500; MCL-07; Active.	3 / 5
	 DB15 (15 pin-macho)	ECG 12; ECG 12S.	10
Esaote	 DB15 (15 pin-macho)	P80; P80 I; P120.	10
Fund. Adib Jatene	 AAMI (6 pin - Macho)	Todos os modelos.	3 / 5

Fukuda	 DB15 (15 pin-macho)	FJC-7110; FCP-2155; FX-2111; FX-4010.	10
Funbec	 DIN 60°	ECG 3.	3 / 5
HW (Heart Ware)	 DB15 (15 pin-macho)	ECG 98; Ergo; ErgoMet.	10
Hill Med	 DB15 (15 pin-macho)	HM-Specialist.	10
HP / Philips / Agilent	 AAMI (6 pin - Macho)	Viridia A1; Viridia A2; 43100A; 43110A; 43120A; Preamp 350;	3 / 5
	 DB15 (15 pin-macho)	Pagewriter (novas versões); XLI; XLS; M1700; M1701; M1770; M1771; M1772; Pagewriter 100; Pagewriter 200; Pagewriter 200i; 300PI.	10
	 HP (12 pin - Macho)	C3; Todos da série M; Intellivue.	3 / 5
Indumed / Mindray	 AAMI (6 pin - Macho)	Vital Line; Vital Wave; Vital Lite; MEC 1000; PM 8000; PM 9000.	3 / 5
Instramed	 SV500	Miniscope.	3 / 5
	 AMP (8 pin - Fêmea)	Miniscope I; Miniscope II.	3 / 5
	 AAMI (6 pin - Macho)	Inmax.	3 / 5
	 SV800	HS 03.	3 / 5
	 AMP (8 pin - Fêmea)	HS 02.	3 / 5
Invivo Research	 AAMI (6 pin - Macho)	Millenia; Millennia; Omni Track; M6; M12.	3 / 5
Marquette Hellige / GE Medical Systems	 AAMI (6 pin - Macho)	Desf. 900; 1200; 3000A; 1500.	3 / 5

	 DB15 (15 pin-macho)	MicroSmart; MAC 500; MAC 100; MAC 1200; MAC 1200 ST; EK36; EK41; EK43; EK51; EK53; EK403; EK413; EK512; EK512P; EK512E; CardioSmart; CardioSys; PPG EK51; PPG EK56; 12-Lead Multilink cable for GE Eagle; Solar; Dash monitors; Cath Lab Systems.	10
Medical Data Eletronics (MDE)	 AAMI (6 pin - Macho)	Escort 100; Escort 101; Escort 102; Escort 200; Escort Prism.	3 / 5
MRL	 AAMI (6 pin - Macho)	360SL/450SL; Porta Pack 80; 90/1000 – Isolated.	3 / 5
Medicor	 DB15 (15 pin-macho)	ER – 65.	10
Medison	 DB15 (15 pin-macho)	EKG 2000.	10
Moriya	 AAMI (6 pin - Macho)	M3000; M9100; M2000.	3 / 5
Mortara	 DB15 (15 pin-macho)	ELI 100; ELI 200; ELI 50.	10
Nellcor / Mallinckrodt	 SV500	NPB–4000; NPB–3900; N–1000; N-200 CLock; N250-C-Lock; 791018.	3
Nihon Kohden	 DB15 (15 pin-macho)	CardioFAX ECG 6511; ECG 6511S; ECG6551; ECG6951D; ECG6851C; ECG8110P.	10
		Cardiofax 5151A; Cardiofax 5403A; Cardiofax 5503A; Cardiofax 5504A; Cardiofax 6151A; Cardiofax 6306A; Cardiofax 6353A; Cardiofax 6543A.	10
	 AAMI (6 pin - Macho)	BSM-7102A; BSM-8301A; BSM-8302A; BSM-8500A; BSM-8502A; BSM-8800A; Life Scope 12; Life Scope 6; Life Scope 7; Life Scope 8; Life Scope 9; MU-831RA; OEC-6105A; OEC-6102A; OEC-7102A; OEC-8108A; TEC-5200A; OMP-7201C; OMP7203C.	3 / 5
Nonicard	 SV500	FAJ.	5
Novamatrix	 AAMI (6 pin - Macho)	902; 903; 910; Life Defense 1.	3 / 5
Schiller	 DB15 (15 pin-macho)	AT1; AT2; AT3; AT4; AT5; AT6; AT10; AT60; CS100; CS200; CS6; AT104; AT101; AT102.	10
Physio Control / Medtronic	 AAMI (6 pin - Macho)	LifePak 5; LifePak 6(S); LifePak 7; LifePak 8; LifePak 9; LifePak 10; Physio 1440/R; VSM-1 (ESF); VSM-2 (ESF); VSM3.	3 / 5
Siemens Eletromedical	 DPK1	1200 Series; 400 Series; 500 Series; 700 Series; 900 Series; XL7000; XL9000.	3 / 5

	 AAMI (6 pin - Macho)	Sirecust DU; Sirecust 630; Sirecust 610; Sirecust 300DU.	3 / 5
	 DB15 (15 pin - fêmea)	Cardiostat 1 / 3 T / 703 / 710; Cardina C; Cardirex 3T/4/6T/62/701/703; Cathcor; Megacard Elema A/B / R / R/E; Mingodoc; Mingocard 3 / 4 / 7; Mingograph 4 / 62 / 81; Mingograph 82 / 740; Mingolog 7; Mingophon 3/7/310/410/420/710/720; Miracard 1; Sirecard 440/440S/460/460S/740; Siereg; Sirecard 803; Sirecust 701.	10
Takaoka	 AAMI (6 pin - Macho)	9100; 9200; Cardiotac; Oxifast.	3 / 5
TEB	 DIN 60°	M10; M10B; Med1.	3
		M10A; M10C.	5
	 DB15 (15 pin-macho)	M10C; C10; ECGPC; APEX 1000.	10
Telemed	 DB15 (15 pin-macho)	Health Tech.	10
Welch Allyn	 AAMI (6 pin - Macho)	Atlas; Encore; Propaq 100 Series; Propaq 200 Series; Propaq CS; Protocol; UltraSmart.	3 / 5
Zool	 AAMI (6 pin - Macho)	D900; PD1200; PD1400; M-Series.	3 / 5